**DUAL LUMEN CATHETER\*  
INSTRUCTIONS FOR USE  
(Peel Away Guidewire Introducer)  
FOR SINGLE USE ONLY**

**CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

**PRODUCT DESCRIPTION OF THE SYSTEM**

The Dual Lumen Catheter is designed to allow for the insertion of two guidewires into a single vessel, preventing them from crossing or entangling, and then allowing them to separately advance into different vessels. The Dual Lumen Catheter consists of two lumens: one lumen allows the guidewire to advance along its length, and the second lumen allows separation and release of the second guidewire via two skives located 35 cm apart and peeling along its length. A radiopaque marker is located at the proximal skive to aid in positioning of the catheter. The Dual Lumen Catheter has a working length of 90 cm and is compatible with a 9 FR introducer sheath. The proximal end of the catheter consists of a female luer connector.

**INTENDED USE**

The Dual Lumen Catheter is intended for use during a two-guidewire procedure.

**CONTRAINDICATIONS**

The Dual Lumen Catheter is contraindicated in patients in whom radiographic contrast or imaging agents are contraindicated.

**POTENTIAL COMPLICATIONS**

Potential complications with this procedure may include, but are not limited to:

- |                           |                   |
|---------------------------|-------------------|
| Vessel perforation        | Vessel dissection |
| Hematoma at site of entry | Wound infection   |

**WARNINGS AND PRECAUTIONS**

- The Dual Lumen Catheter shall be used only by physicians trained in vascular surgery, interventional radiology, or interventional cardiology.
- The Dual Lumen Catheter should be handled with care to avoid excessive bending or kinking of the catheter, which may result in an inability to properly advance the guidewires.
- To avoid insertion difficulties or tissue injury, never advance this device without the use of fluoroscopic guidance.
- The Dual Lumen Catheter is not intended to be used to infuse radiographic contrast agents.
- The Dual Lumen Catheter is designed for single use only. Do not reuse or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

***READ ALL INSTRUCTIONS THOROUGHLY  
BEFORE USE.***

\* US Patent: 6,440,161

Other US and Foreign Patents Pending

**INSTRUCTIONS FOR USE**

1. Perform standard procedure to access vessel.
2. Flush lumen.
3. Advance a 0.89 mm (0.035") guidewire to the desired location for the Dual Lumen Catheter placement.
4. Advance the Dual Lumen Catheter over this guidewire via the inner reinforced catheter lumen.
5. Advance the Dual Lumen Catheter so that the radiopaque marker (located at the proximal skive of the outer catheter lumen), is in the desired location.

6. Remove the guidewire used for Dual Lumen Catheter placement.
7. Insert the first guidewire through the inner reinforced catheter lumen.
8. Insert the second guidewire into the distal skive until it exits the proximal skive. This will keep the two guidewires separated.
9. To remove the second guidewire, retract the inner reinforced catheter lumen, and the guidewire will separate as the outer catheter lumen peels away.

**WARNING: CATHETER ADVANCEMENT SHOULD BE PERFORMED UNDER FLUOROSCOPIC GUIDANCE. DO NOT USE EXCESSIVE FORCE TO ADVANCE OR WITHDRAW THE CATHETER WHEN RESISTANCE IS ENCOUNTERED.**

**THE DUAL LUMEN CATHETER IS DESIGNED FOR SINGLE USE ONLY.**

**STORAGE**

Store the Dual Lumen Catheter in a cool, dry place.

**STERILIZATION**

The Dual Lumen Catheter is a single use item and shipped sterile. Contents sterile and non-pyrogenic only if the package is not opened, damaged or broken.

**WARNING: The Dual Lumen Catheter is designed for single use only. Do not reuse or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.**

**INSTRUCTIONS FOR RETURNING PRODUCT**

In the event the device must be returned for any reason, place the Dual Lumen Catheter into its original package and shipping box. Contact Customer Service inside the US at 800-983-2284 or Outside US at +1 (949) 595-7200 to receive a return goods authorization number (RGA), and ship the device to the address listed on the label.

**SYMBOL LEGEND**

SYMBOL	DESCRIPTION
	Product expiration date
	Batch Code, Lot number or work order number for the product
	Method of sterility is Ethylene Oxide
	Caution. Consult accompanying documents.
	Device is intended for use one time only. Do not re-use or re-sterilize.
	CE Mark symbol with Notified Body registration number (NSAI)
	Catalog Number
	Consult Instructions for Use.
	Authorized Representative in the European Community (NSAI)
	Manufacturer
	Do not use if package is damaged.
	Keep Dry.

SYMBOL	DESCRIPTION
	Keep Away From Heat.

## **MANUFACTURER ADDRESS**

Endologix, Inc.  
 11 Studebaker  
 Irvine, CA 92618-2013 USA  
 Tel: Inside US (800) 983-2284  
 Outside US (949) 595-7200  
 Fax: (949) 830-4463

**EMERGO EUROPE**  
 Molenstraat 15  
 2513 BH, The Hague  
 The Netherlands  
 Tel: +31 (0)70 345 8570  
 Fax: +31 (0)70 346 7299

©2010 Copyright. All Rights Reserved.  
 Part Number C00078 Revision F

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©Copyright 2010. Tous droits réservés.
©2010 Copyright. Tutti i diritti riservati.
©2010 Copyright. Reservados todos los derechos.
©2010 Copyright. Alle Rechte vorbehalten.
©2010 Copyright. Todos os direitos reservados.
©2010 Copyright. Alle rettigheder forbeholdes.